

Product datasheet for MR217835

Galr3 (NM_015738) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Galr3 (NM_015738) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Galr3
Synonyms:	Galnr3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR217835 representing NM_015738 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGACATCCAGAACATTTTCGCTGGACAGCCAGGGAGCGTAGGGCTGTGGCAGTGCCGGTGGTCT
TTGCCCTCATCTTCTGTTGGGCATGGTGGTAATGGGCTGGTGTGGCCGTGCTTCTGCAGCCTGGCCC
AAGTGCCCTGGCAGGAACCGGGCAGTACCACGGATCTTTCATCCTCAACTGGCCGTGGCCGACCTTTGC
TTCATCCTGTGCTGCGTGCCCTCCAGGCGGCCATCTACACACTGGATGCCTGGCTTTTGGGGCTTTCG
TGTGCAAGACGGTACATCTGCTCATCTACCTACCATGTATGCCAGCAGCTTACCCTGGCGGCCGTCTC
AGTGGATAGGTACCTGGCGGTGCGACACCCGCTGCGCTCCCGGGCCCTCCGCACCCCGCGCAACGCGCGC
GCCGCCGTGGGGCTCGTGTGGCTGCTGGCGGCGCTCTTTTCCGCGCCCTACCTAAGCTACTACGGCACGG
TGCGCTACGGCGCGCTCGAGCTCTGCGTGCCCGCTGGGAGGACGCGCGGCGCGCGCTGGACGTGGC
CACCTTCGCCGACGGCTACCTGCTGCCGGTGACCGTGGTGGAGCTGGCCTACGGACGCACGCTGTGCTTC
CTGTGGGCCCGCGTGGGTCCCGCGGGTGGCGGGCGGCAGAGGCACGCAGACGGGCAACCGGCCGCGCG
GGCGGCCATGCTGACGGTGGCCGCGCTACGCACCTGCTGGGGCCCGCACACGCGCTCATCTCTG
CTTCTGGTACGGCCGCTTCGCCCTCAGCCCGGCCACCTACGCCTGCTGCTGGCTTCGCACTGCCTCGCC
TACGCCAACTCCTGCCTCAACCCGCTCGTCTATTCCCTTGCCCTCCGCCACTTCCGCGCGGCTTCCGCC
GCCTGTGGCCCTGCGGCCACCGCCGCCACCGCCACCACCACCACCGCCTTCATCGAGCCCTCCGTGCTGT
CCAGCCGCGCTTTCGGGCCCGCGGCTATCCCGGCGACGCCAGACCTCGCGGTTGGAGTATGGAGCCC
AGAGGGGATGCCTTGCGCGGTGGAGAGACTAGACTAACCTGTCCGCCAGGGGACCGCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



[View online »](#)

Protein Sequence: >MR217835 representing NM_015738
Red=Cloning site Green=Tags(s)

MADIQNISLDSPGSVGAVAVPVVFALIFLLGMVGNGLVLAVLLQPGPSAWQEPGSTTDLFILNLAVADLC
 FILCCVPFQAAIYTLDAWLFGAFVCKTVHLLIYLTMYASSFTLAAVSVDRYLAVRHPLRSRALRTPRNAR
 AAVGLVWLLAALFSAPYLSYYGTVRYGALELCVPAWEDARRRALDVATFAAGYLLPVTVVSLAYGRTLCF
 LWAAVGPAGAAAAEARRRATGRAGRAML TVAAL YALCWGPHHALILCFWYGRFAFSPATYACRLASHCLA
 YANSCLNPLVYSLASRHFRRFRRLWPCGHRHRHHHRLHRALRRVQPASSGPAGYPGDARPRGWSMEP
 RGDALRGGETRLTLSARGPQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9019_d01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_015738


ORF Size: 1110 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

- OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
- RefSeq:** [NM_015738.2](#), [NP_056553.2](#)
- RefSeq Size:** 1245 bp
- RefSeq ORF:** 1113 bp
- Locus ID:** 14429
- UniProt ID:** [O88853](#)
- Cytogenetics:** 15 37.7 cM
- MW:** 40.9 kDa
- Gene Summary:** Receptor for the hormone galanin and spexin-1.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR217835