

Product datasheet for **MR217917L4V**

Armcx2 (NM_026139) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Armcx2 (NM_026139) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Armcx2
Synonyms:	3230401N03Rik; AI043003; ALEX2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026139
ORF Size:	2352 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217917).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026139.4 , NP_080415.3
RefSeq Size:	3790 bp
RefSeq ORF:	2355 bp
Locus ID:	67416
UniProt ID:	Q6A058
Cytogenetics:	X E3
Gene Summary:	May regulate the dynamics and distribution of mitochondria in neural cells.[UniProtKB/Swiss-Prot Function]



[View online »](#)