

Product datasheet for **MR217920L4V**

Csf2ra (NM_009970) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Csf2ra (NM_009970) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Csf2ra
Synonyms:	CD116; Csfgmra; GM-CSF-Ra; GM-CSFR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009970
ORF Size:	1167 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR217920).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009970.2 , NP_034100.2
RefSeq Size:	1850 bp
RefSeq ORF:	1167 bp
Locus ID:	12982
UniProt ID:	Q00941
Cytogenetics:	19 56.89 cM
Gene Summary:	Low affinity receptor for granulocyte-macrophage colony-stimulating factor. Transduces a signal that results in the proliferation, differentiation, and functional activation of hematopoietic cells.[UniProtKB/Swiss-Prot Function]



[View online »](#)