

Product datasheet for **MR218021L3V**

Rnf149 (NM_001033135) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rnf149 (NM_001033135) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rnf149
Synonyms:	1600023E10Rik; AV037261; C76763; Gm15832; Greul4; OTTMUSG00000026591
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001033135
ORF Size:	1182 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218021).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001033135.3 , NP_001028307.2
RefSeq Size:	2440 bp
RefSeq ORF:	1185 bp
Locus ID:	67702
UniProt ID:	Q3U2C5
Cytogenetics:	1 B
Gene Summary:	E3 ubiquitin-protein ligase. Ubiquitinates BRAF, inducing its proteasomal degradation. [UniProtKB/Swiss-Prot Function]



[View online »](#)