

## Product datasheet for **MR218035L3V**

### Eid3 (NM\_025499) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | Eid3 (NM_025499) Mouse Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | Eid3   |
| Synonyms:                 | 1700027M21Rik; NS4EB   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_025499  |
| ORF Size:                 | 930 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(MR218035).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_025499.2</a> , <a href="#">NP_079775.2</a>  |
| RefSeq Size:              | 1305 bp  |
| RefSeq ORF:               | 1128 bp  |
| Locus ID:                 | 66341  |
| UniProt ID:               | <a href="#">Q3V124</a>   |
| Cytogenetics:             | 10 C1  |



[View online »](#)

**Gene Summary:**

Tissue-specific component of the SMC5-SMC6 complex, a complex involved in repair of DNA double-strand breaks by homologous recombination. The complex may promote sister chromatid homologous recombination by recruiting the SMC1-SMC3 cohesin complex to double-strand breaks. The complex is required for telomere maintenance via recombination and mediates sumoylation of shelterin complex (telosome) components (By similarity). [UniProtKB/Swiss-Prot Function]