

Product datasheet for **MR218072L4V**

C1rl (NM_181344) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C1rl (NM_181344) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	C1rl
Synonyms:	C1r-LP; C1rl1; C1rlp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_181344
ORF Size:	1446 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218072).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_181344.5 , NP_851989.3
RefSeq Size:	3010 bp
RefSeq ORF:	1449 bp
Locus ID:	232371
UniProt ID:	Q3UZ09
Cytogenetics:	6 F2
Gene Summary:	Mediates the proteolytic cleavage of HP/haptoglobin in the endoplasmic reticulum. [UniProtKB/Swiss-Prot Function]



[View online »](#)