

## Product datasheet for **MR218079**

### Zfp467 (NM\_001085416) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Zfp467 (NM\_001085416) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Zfp467  
**Synonyms:** 1190001I08Rik; AB041616; AI426952; AV017623; EZI; MNCb-3350; Znf467  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR218079 representing NM\_001085416  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGAGAGACCCTGGAGGCCCTGAACCTCTGGGATTCTCAGTGGGACAGCCAGAGATGGCTCCCCAGA  
GTGAGCCAGGGATGGGTTTCAGTAATGCCAGGAGAAGATGTCATCCAGAGGAGAGACACTGCCTC  
CTGCTCAGGGCATGAGACTCCTGGCCAGAAAGAAGGCATCCACACAGAACAAGCTGAAGCTCCTGTATG  
GGCAGCCAGGCCAGTACCCACAGAAGGCCGAACCTGCAGGCTCTGTCCCAGTCTTTCCCACTCTGA  
CTTGGGTCAGTTGGCTGTACGTACAACTGGAGCCAGGGACTCCGGGAGCAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR218079 representing NM\_001085416  
Red=Cloning site Green=Tags(s)  
MRETLEALNSLGFVSGQPEMAPQSEPRDGFSAQEKMSRGESTLHSCSGHETPGQKEGIHQAEAPCM  
GSQASTPQKAEPAGSVPVLFPLLTVVSWLLRTNWSQGLREH

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9042\\_a02.zip](https://cdn.origene.com/chromatograms/mm9042_a02.zip)

**Restriction Sites:** SgfI-MluI



[View online >](#)

**Cloning Scheme:**


**ACCN:** NM\_001085416

**ORF Size:** 333 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001085416.1](#), [NP\\_001078885.1](#)

**RefSeq Size:** 3179 bp

**RefSeq ORF:** 336 bp

**Locus ID:** 68910


UniProt ID: [Q8JZL0](#)

Cytogenetics: 6 B2.3

MW: 12.5 kDa

**Gene Summary:** Transcription factor that promotes adipocyte differentiation and suppresses osteoblast differentiation in the bone marrow. Enhances the osteoclast-supporting ability of stromal cells. Binds with STAT3 the consensus sequence 5'-CTTCTGGGAAGA-3' of the acute phase response element (APRE). Transactivates several promoters including FOS, OSM and PPARG. Recruits a histone deacetylase complex.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR218079