

Product datasheet for **MR218106L3V**

Pde4a (NM_019798) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pde4a (NM_019798) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pde4a
Synonyms:	D9Ertd60e; Dpde2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019798
ORF Size:	1830 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218106).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019798.5 , NP_062772.3
RefSeq Size:	4150 bp
RefSeq ORF:	1833 bp
Locus ID:	18577
UniProt ID:	O89084
Cytogenetics:	9 7.73 cM
Gene Summary:	Hydrolyzes the second messenger cAMP, which is a key regulator of many important physiological processes.[UniProtKB/Swiss-Prot Function]



[View online »](#)