

## Product datasheet for **MR218126L3V**

### **Ehbp111 (NM\_053252) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ehbp111 (NM_053252) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ehbp111
Synonyms:	G430002G23Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053252
ORF Size:	2631 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218126).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053252.3</a> , <a href="#">NP_444482.2</a>
RefSeq Size:	3888 bp
RefSeq ORF:	2634 bp
Locus ID:	114601
UniProt ID:	<a href="#">Q99MS7</a>
Cytogenetics:	19 A
Gene Summary:	May act as Rab effector protein and play a role in vesicle trafficking (By similarity). Involved in apical-directed transport in polarized epithelial cells; the functions seems to implicate Rab8, BIN1 and possibly DNM1 (PubMed:26833786).[UniProtKB/Swiss-Prot Function]



[View online »](#)