

## Product datasheet for **MR218176L3V**

### **Ss18 (NM\_001161371) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ss18 (NM_001161371) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ss18
Synonyms:	D130059H17; Ssxt; Syt
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001161371
ORF Size:	939 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218176).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001161371.1</a> , <a href="#">NP_001154843.1</a>
RefSeq Size:	2795 bp
RefSeq ORF:	942 bp
Locus ID:	268996
Cytogenetics:	18 A1
Gene Summary:	Appears to function synergistically with RBM14 as a transcriptional coactivator. Component of SWI/SNF chromatin remodeling subcomplex GBAF that carries out key enzymatic activities, changing chromatin structure by altering DNA-histone contacts within a nucleosome in an ATP-dependent manner.[UniProtKB/Swiss-Prot Function]



[View online »](#)