

## Product datasheet for **MR218183L4V**

### Ccl8 (NM\_021443) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ccl8 (NM_021443) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ccl8
Synonyms:	1810063B20Rik; AB023418; HC14; MCP-2; Mcp2; Scya8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021443
ORF Size:	291 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218183).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021443.3</a> , <a href="#">NP_067418.1</a>
RefSeq Size:	517 bp
RefSeq ORF:	294 bp
Locus ID:	20307
UniProt ID:	<a href="#">Q9Z121</a>
Cytogenetics:	11 49.91 cM
Gene Summary:	Chemotactic factor that attracts monocytes. This protein can bind heparin (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)