

Product datasheet for **MR218190L3V**

Tcf25 (NM_025804) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tcf25 (NM_025804) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tcf25
Synonyms:	1100001J13Rik; 1810041K11Rik; D8Ertd325e; mKIAA1049; Nulp1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025804
ORF Size:	1875 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218190).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025804.2
RefSeq Size:	2695 bp
RefSeq ORF:	1878 bp
Locus ID:	66855
UniProt ID:	Q8R3L2
Cytogenetics:	8 72.1 cM
Gene Summary:	May play a role in cell death control. Acts as a transcriptional repressor. Has been shown to repress transcription of SRF in vitro and so may play a role in heart development (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)