

## Product datasheet for **MR218202L4V**

### Taf11 (NM\_026836) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Taf11 (NM_026836) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Taf11
Synonyms:	28kDa; 1110038O14Rik; AI115401; AL024207; AW552398; PRO214; TAF2I; TAFII28
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026836
ORF Size:	636 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218202).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026836.2</a>
RefSeq Size:	1220 bp
RefSeq ORF:	636 bp
Locus ID:	68776
UniProt ID:	<a href="#">Q99JX1</a>
Cytogenetics:	17 A3.3
Gene Summary:	Core TAFII present in both of the previously described TFIID species which either lack or contain TAFII30 (TFIID alpha and TFIID beta respectively).[UniProtKB/Swiss-Prot Function]



[View online »](#)