

Product datasheet for MR218222

Pigp (NM_001159619) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pigp (NM_001159619) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Pigp
Synonyms: Dcrc; Dscr5
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR218222 representing NM_001159619
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAGATTATGTAACAGACCATTGGTGAAGGGGGCGGGGGGAGATGTCTGACGACGATTGTGTTCCCTT
TCCACAGATTGTCTAAAGCCACAGGGAAGATGGTGGAAAATTCACCGTCGCCATTGCCAGAAAGAGCGAT
TTATGGCTTTGTCTTTTCTAAGCTCACAATTTGGCTTCATACTTTATCTTGTGTGGGCTTTCGTTCCC
GAATCTTGTTAAACTCCTTAGGCTTAAGTACTGGCCTCAAAAATACTGGGCGGTTGCATTGCCCGTCT
ACCTCCTTATCACAGTTGTAATTGGTTATGTGCTCTTATTTGGAATCAACATGATGAGTACCTCTCCACT
TGACTCCATTACATACGATCACAGATAACTATGCTAAAAACCAACAGAGGAAAACTACCAAGAGGATGCC
ATCCCAGCATTGAGAGATGTCCCATTAGTGAAGTAAACAAAATGTTTTCTTAGGAGCCAAAGAACTCA
ATACC

ACGCGTACGCGGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR218222 representing NM_001159619
Red=Cloning site Green=Tags(s)

MRLCNRPLVKGGRCLTTIVFLFHRLSKATGKMVENSPLPERAIYGFVFLSSQFGFILYLWVAFVP
ESWLNLSGLTYWPQKYWAVALPVYLLITVVIGYVLLFGINMMSTSPLDSIHTITDNYAKNQQRKNYQEDA
IPALRDVPISEVNKMFFLGAKELNT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001159619

ORF Size: 495 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001159619.1</u> , <u>NP_001153091.1</u>
RefSeq Size:	805 bp
RefSeq ORF:	498 bp
Locus ID:	56176
UniProt ID:	<u>Q9JHG1</u>
Cytogenetics:	16 C4
MW:	19.2 kDa
Gene Summary:	Part of the complex catalyzing the transfer of N-acetylglucosamine from UDP-N-acetylglucosamine to phosphatidylinositol, the first step of GPI biosynthesis. [UniProtKB/Swiss-Prot Function]