

Product datasheet for MR218226

Rpl6 (NM_011290) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rpl6 (NM_011290) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rpl6
Synonyms:	Taxreb107
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR218226 representing NM_011290 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGGTGAAAAGCGCCTGATACAAAGGAGAAGAAACCTGCGGCCAAGAAGGCTGGCAGTGATGCTG
CTGCCTCCCGGCCCGGGCCGCAAGTAGCCAAAAGGTTACCCCTAAGGGTAAAAAGCCCAAGAAGGC
GAAGCCCCATTGCAGCCGAAACCCTGTCTGGTGAGAGGAATTGGCAGGTAATCCCGATCTGCTATGTAC
TCCAGGAAGGCCTGTACAAAAGGAAATACTCGGCTGCCAAGACCAAGGTTGAGAAGAAGAAGAAAAAG
AGAAGGTCCTTGCTACTGTACCCAAAACAGTTGGTGGGACAAGAACGGTGGCACCCGGGTGGTGAAGCT
TCGAAAAATGCCTAGGATTATCCACCGAAGATGTTCCCTCGGAAAGCTGCTGAGCCACGGCAAGAAGCCC
TTCAGCCAGCACGTGAGAAGGCTACGCTCCAGCATCACTCCCGAACTGTCTGATCATCTCACTGGGC
GCCACAGGGGCAAGAGAGTGGTTTTCTGAAGCAGCTGGACAGTGGCTTGTCTGTGACTGGGCCTCT
TGTATCAACAGAGTTCTCTGCGCAGAACACACCAGAAGTTTGTATTGCCACCTTACAAAAGTTGAT
ATCAGCGATGTTAAAATCCCAAACACCTGACTGACGCTTACTTCAAGAAGAAGCAGCTGCGCAAGCCCA
GGCATCAGGAGGGCGAGATCTTCGACACAGAGAAGGAGAAATACGAGATTACAGAGCAGCGAAAGGCTGA
CCAGAAAGCTGTGGATTTGCAGATTTTGCCAAAGATTAAGCTGTTCTCAGCTCCAGGGCTACCTGGC
TCTCAGTTTTCCCTGACAAACGGGATGTATCCTCACAACCTGGTCTTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR218226 representing NM_011290
 Red=Cloning site Green=Tags(s)

MAGEKAPDTKEKKPAAKKAGSDAAASRPRAAKVAKKVHPKGGKPKKAKPHCSRNPVLRGIGRYSRSAMY
 SRKALYKRKYSAAKTKVEKKKKKEKVLATVTKTVGGDKNGGTRVVKLRKMPRYPTEDVPRKLLSHGKKP
 FSQHVRRLRSSITPGTVLIILTGRHRGKRVVFLKQLDSGLLLVTGPLVINRVPLRRTHQKFVIATSTKVD
 ISDVKIPKHLTDAYFKKKQLRKPRHQEGEIFDTEKEYEITEQRKADQKAVDLQILPKIKAVPQLQGYLR
 SQFSLTNGMYPHKLVF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_011290

ORF Size: 888 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_011290.5](#), [NP_035420.2](#)

RefSeq Size: 1282 bp

RefSeq ORF: 891 bp

Locus ID: 19988


UniProt ID: [P47911](#)

Cytogenetics: 5 F

MW: 34 kDa

Gene Summary: Component of the large ribosomal subunit.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR218226