

Product datasheet for MR218240L3

Msmb (NM_020597) Mouse Tagged Lenti ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Msmb (NM_020597) Mouse Tagged Lenti ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | Msmb |
| Synonyms: | beta-MSP; PIP; PSP94 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| E. coli Selection: | Chloramphenicol (34 ug/mL) |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR218240). |
| Restriction Sites: | SgfI-MluI |
| Cloning Scheme: | |

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

| | |
|-----------|-----------|
| ACCN: | NM_020597 |
| ORF Size: | 339 bp |



[View online »](#)

