

Product datasheet for **MR218253**

Lrrc63 (NM_027581) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lrrc63 (NM_027581) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Lrrc63
Synonyms:	4921509B22Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR218253 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAAAAGAGCATTATCGCTCACTCAAAACCCAGCAGCATCCACTACTTCTGCGAAGACCCTTACCTC
 CAAAGCTCCCTAAAGTACCTTTATTAATAAAGAAAGTTTCTGATAGTACATACAGTTAAGGCTAGAGAACC
 GAAGGAGGAGACCTATGGAAAGTTGATCCTGATGAAAGTGCTTTCGTGAAAATGGATCAAAATGTTTCT
 CAATCACAAGACCTACCAGCACTCAAACCATCTGTCTTGACTATGACCAATAAGGCACTCGAAAGAGTGG
 TTAACATTTTTCCACATCAACACCGTGTTCGCCGTTTTCACTGGAAGGTCCCAAAGACACCGCTGGTTC
 CATGTTCAATCCAAGATGTCTGTCGGCGTCTCCAGAATTTTTAAGAATACGCTATCCCAATGAAAAA
 AGAGCCAAGAAGTCTCGAAAAAGGAAGTCAAAGAAGATGGGTATCTCACTAAAAAGACATTTAACACTT
 TAGTTCCTTCCAAGGAGTTTTCAAACCTTCACCCACTTCTACTCTAGATTCATCGCATCTAAGTTTAA
 GCCACTAGGTACCCAGTTTCACCTGCACCAAAGATTTCAAGCTATAGTCAAGAGCTCCCTCTTCTCAGA
 GATGTAAGAGGGCAGCACTCAGCCCTGTACCTAGCACTTCATCTGTAATACCTGAACCGCAGTGGTCTG
 AACGGACACACCTTCGGCAGCCCAAGTAAAGTGGTACTGAACACTGGCTCACTCCCGCCAGCCCGTTC
 TCTGCCAACGCCTGTTTTGCCAAGAAAACCTCTAAGACAGGCTATGATAGAAAATGCTGCTGCTGCCGCT
 GCC
 CCACCACCGCTGCAGTGAAGTGGAAAAACAGAAGCCACAAACCGCTGAAACTGTTGAGAGAACCACAG
 ACAAATAGGCCAGATGCCCATGTCTTACGAGGGGAAGGCTTTAAGGCTGTTGCAGCAACACGCTCTGAA
 ACAGTACTGGCGCTGACCACCTGGCCATAATAAATTGCCAAATATATGGAAGAAATGCACTTAATCTTA
 AGGGATTCTTTCTAGCGAACTGCCAGACTTAACACCTGTGGCTTTTCAATTGGTATACCTTAATTTATC
 ATTTAATGACTTAAACCAATTTCCATAGAAATACTTTATCTTCAAACCTACAAGTACTGAAATTGAGA
 AATAACCTATCAAAGAAATTCCTTCTGAAATCCATCTATTGACATACCTTAGGATTTTCTCCATTGCCT
 TTAATTATATTACAAAGCTTCTGATGGGTTATTTGCCTGAATTATCTTGAGGAGCTCGATGTTTCCTA
 CAATGAAATTGAAAATATTTCAAATGAAATCCAGAAACTCAGATCACTTGAAAAACTGATTGTTGATGGG
 AATCCCATAACTCTTTCCACCTGGGATTTTGAAGCTTAACTTGATCAAATTACAAATTGAAAATACAT
 TCACTTGTCTCAGTTTTGGTTAGAGAGTTCCTGGAATGACCCCAACAACCTACTCAAATTTGCTCCTT
 GTTCTTGTAAAAATAAACTGTTGGATTACATCCAGATGCAGTTCGGAAGTCTCTTAAAAGCACATCT
 GAATGTGATTGGTGCCATGGTCCAAGTTTGGTGAAGTTTCCGCATCATCCGATCCTGCGATATCTTTG
 GAGTATACATGTCCCATTATGTTCCATGTCTGCTCCTACTTGTATCGGGACATCAGGGAACCAG
 CTTTGTTTTAGAAGGTTTCCCTAGTAGAAGAATAGCTCTACATATGGACTGGGTAAAAGAGAGCAAAGTC
 AGCAATGTCTCGTTCTATTTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR218253 protein sequence
 Red=Cloning site Green=Tags(s)


MKKSIIYRSLKTQQHPLLRPLPPKLPKVPLLLKKKVRIVHTVKAREPKEETYGKFDPEESAFVKMDQNVN
 QSQRPTSTQICLDYDHKALERVVNIFFPHQHRVRRFHWKVPKSTAGSMFIPRCLSASSRIFKNTLSQMKK
 RAKKSSKKEVKEDGYLTKKTFNTLVLSKEFSKPSPTSYSRFIASKFKPLGTQFHPAPKISGYSQELPLLR
 DVKRAALSPVPSTSSVIPEPQWSEKTHPSAAPSQVVLNTGSLPPARSLPTVLPKPLRQAMIENAAAAA
 AAAAAAAAAAAAAATTTTTTAAATTAAVSGKTEAHKPPETVQRTTRTIGPDAHVLRGEGFKAVAAATRSE
 TVLALTTLAIINCQIYGRNALNLKGFFLANCPDLTPVAFQLVYLNLNFNDLNQFPIEILYLQNLQVLKLR
 NNPIKEIPSEIHLTYLRIFSIAFNYITKLPDGLFCLNYLEELDVSNEIENISNEIQKLRSLKLVLDG
 NPITSFPFGILKLNLIKQIENFTCSQFWLESSWNPQQLTQICSLFLVKNKLLDYPDAVRKSLKSTS
 ECDWCHGPKFGEFRIIRSCDIFGVSHVPIMFHVCSSTCYRDIRETSFVLEGFPSRRIALHMDWVKESKV
 SNVSFYL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_027581

ORF Size: 1914 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_027581.1](#), [NP_081857.1](#)

RefSeq Size: 2519 bp

RefSeq ORF: 1917 bp


Locus ID: 70859

UniProt ID: [A6H694](#)

Cytogenetics: 14

MW: 71.5 kDa

Product images:



Circular map for MR218253