

Product datasheet for **MR218269L4V**

Lrrc10 (NM_146242) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lrrc10 (NM_146242) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lrrc10
Synonyms:	D330003I11Rik; Hrlrrp; Sardin1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_146242
ORF Size:	825 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218269).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_146242.2 , NP_666354.1
RefSeq Size:	1428 bp
RefSeq ORF:	825 bp
Locus ID:	237560
UniProt ID:	Q8K3W2
Cytogenetics:	10 D2
Gene Summary:	May play important roles in cardiac development and/or cardiac function.[UniProtKB/Swiss-Prot Function]


[View online »](#)