

Product datasheet for **MR218277L3V**

Baalc (NM_080640) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Baalc (NM_080640) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Baalc
Synonyms:	2810457D07Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_080640
ORF Size:	435 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218277).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080640.5
RefSeq Size:	4177 bp
RefSeq ORF:	438 bp
Locus ID:	118452
UniProt ID:	Q8VHV1
Cytogenetics:	15 15.18 cM
Gene Summary:	May play a synaptic role at the postsynaptic lipid rafts possibly through interaction with CAMK2A.[UniProtKB/Swiss-Prot Function]



[View online »](#)