

Product datasheet for **MR218282L4V**

Als2cl (NM_001146059) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Als2cl (NM_001146059) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Als2cl
Synonyms:	79H19D; 79H19E; BC022645; D930044G19Rik; mRn.49018
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001146059
ORF Size:	2859 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218282).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001146059.1
RefSeq Size:	5106 bp
RefSeq ORF:	2859 bp
Locus ID:	235633
UniProt ID:	Q60I26
Cytogenetics:	9 60.79 cM
Gene Summary:	Acts as a guanine nucleotide exchange factor (GEF) for Rab5 GTPase. Regulates the ALS2-mediated endosome dynamics (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)