

Product datasheet for **MR218283**

Rnpep (NM_001159624) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rnpep (NM_001159624) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rnpep
Synonyms:	MGC29229
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR218283 representing NM_001159624
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGAGCGGCCGACCAGGCAACTACAGCGCGGCCGCCCGCGGCCACTGCACTCGGCCAAGCCGTGG
 ACGTGGCCTCGGCTCCAGCTTCCGCGCCTTCAGATCCTGCACCTGCACCTGGACCTGCGGGCCGAGTT
 CGGGCCTCCGGGCCGGGTCCGGGAAGCCGCGGATTGAGCGGCACGGCGACGCTGGAGCTGCGTTGTTTG
 CTTCCCGAGGGCGCTCGGAGCTGCGGCTCGACTCGCACAGCTGCTTGGAGGTGACGGCGGCGACGCTGC
 GACGGGGACAGCCGGGAGACCAGCAGGCCCGGCAGAGCCGTGCCTTTCCACACACAGCCCTTCTCGCA
 CTATGGCCAGGCCCTGTGCGTGGCCTTCCGGCAGCCCTGCGGTGCCGCCGATCGTTGAGCTGGAGCTC
 ACCTACCGCGTCGGGGAGGGACCCGGGTTTGTGTTGGCTCCTGAGCAGACAGCAGGAAAGAAGAAAC
 CCTTTGTCTATACCCAGGGTCAGGCGGTCTTAACCGAGCCTTTTCCCTTGCTTTGACACTCCAGCAGT
 GAAATGCACGTATTCGGCTCTTATTGAGGTTCCAGATGGCTTCCAGCTGTGATGAGTGCAGACACCTGG
 GAGAAGAGAGGACCAAATAAGTTCTTCCAGATGAGCCATCCCATCCCTCCTATCTGATAGCTTTGG
 CCATTGGAGATCTGGCTCAGCTGAGGTTGGACCCAGGTATGACTTGCTCTTTATGCCACCGTCTTTCC
 ATTTGGAGGAATGGAAAATCCTTGCTGACCTTTGTCACTCCGTGCTGTAGCAGGGGACCGCTCGTTA
 GCCGACGTTATCATCCATGAGATCTCCACAGTTGGTTTGGGAACCTGGTACCAATGCCAACTGGGGGG
 AGTTCTGGCTCAATGAAGGCTTACCATGTATGCCAGAGGAGGATCTCCACCATCTGTTTGGAGCTGC
 TTACACTTGCTTGGAGGCTGCCACAGGGAGGGCTCTGCTGCGTCAGCACATGAATGTCTCCGGAGAGGAA
 AACCCACTCAACAAGCTCCGGGTGAAGATTGAACCAGGTGTGGACCCGATGACACCTACAATGAGACCC
 CCTACGAGAAAGTTACTGCTTTGTCTCATACCTGGCCACCTGGTGGGTGATCAGGATCAGTTTGACAA
 GTTTCTCAAGGCCTATGTAGATGAGTTAAATTCCAAAGTATCTTGGCTGAAGACTTTTGGAAATCTAC
 CTGGAGTATTTCCCTGAGTTGAAGAAGAAGGGAGTAGATTCCATTCCAGTTTTGAGTTTATCGATGGC
 TGAACACCCCTGGATGGCCCCCTACCTCCCGACCTCTCCCTGGGGACTCACTCATGAAGCCAGCTGA
 AGAGCTGGCAGAAGTGTGGGTGACCTCAGAGCCAGACATGCAAGCCATTGAAGCCGTGGCCATCTTACC
 TGGAAAGACCTACCAGCTTGTCTACTTCTAGATAAGATCTTGCAGAAGTCTCCTCTCCCTCTGGGAATG
 TGAAAAAAGTGGGGAGACATACCCAAAGATTTCCAATGCTCAGAATGCTGAGTTACGGCTAAGTGGGG
 CCAATCATCCTAAGAATGACTACCAAGAGGAATCCAGAAAGTGAAGGACTTCTGCAGAGCCAGGGA
 AAGCAGAAGTACTCTGCTCTGTACCACGCCATGATGGGTGGCAGTGAAGTGGCCCGGACCTGGCCA
 AGGACACCTTTGCCGCCACTGCTTCCAGCTCCACAGCAACGTCGTAACATATGCCAGCAGATCTGGC
 GCCAAGGACAGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>MR218283 representing NM_001159624
 Red=Cloning site Green=Tags(s)

MESGGPGNYSAARRPLHSAQAVDVASASSFRAFEILHLHLDLRAEFGPPGPGPSRGLSGTATLELRCL
 LPEGASELRDLSHSCLEVTAAATLRRGQPGDQQAPEVPFHTQPFVSHYQALCVAFRQPCGAADRFELEL
 TYRVGEGPGVCWLAPEQTAGKKKPFVYQQAVALNRAFFPCFDTPAVKCTYSALIEVPDGFATVMSADTW
 EKRGPNKFFQMSHPIPSYLIALLAIGDLASAIEVGPYDILLFMPSPFFGGMENPCLTFVTPCLLAGDRSL
 ADVIIHEISHSWFNLVTNANWGEFNLNEGFTMYAQRRIITILFGAAYTCLAAATGRALLRQHMMNVSGEE
 NPLNKLRVKIEPGVDPDDTYNETPYEKGYCFVSYLAHLVGDQDQDFKFLKAYVDEFKQFQSIILAEDFLEFY
 LEYFPELKKKGVDIPGFEFDRWLNTPGWPPYLPDLSPGDSLMKPAEELAEWVWSEPDMQAIIEVAIAIST
 WKTYQLVYFLDKILQKSPLPPGNVKKLGETYPKISNAQNAELRLRWGQIILKNDYQEEFQKVKDFLQSQG
 KQKYTLPLYHAMGGSEMARTLAKDTFAATASQLHSNVVNYVQQILAPKDS


TRTRPLEQKLISEEDLAANDILDYKDDDDKVV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001159624

ORF Size: 1833 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001159624.1 , NP_001153096.1
RefSeq Size:	2184 bp
RefSeq ORF:	1836 bp
Locus ID:	215615
UniProt ID:	Q8VCT3
Cytogenetics:	1 E4
MW:	68.5 kDa
Gene Summary:	Exopeptidase which selectively removes arginine and/or lysine residues from the N-terminus of several peptide substrates including Arg(0)-Leu-enkephalin, Arg(0)-Met-enkephalin and Arg(-1)-Lys(0)-somatostatin-14. Can hydrolyze leukotriene A4 (LTA-4) into leukotriene B4 (LTB-4) (By similarity).[UniProtKB/Swiss-Prot Function]