

Product datasheet for **MR218296A1V**

Mouse Slirp (NM_026958) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Slirp (NM_026958) AAV Particle
Tag: Myc-DDK
Symbol: Slirp
Synonyms: 1810035L17Rik
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR218296 representing NM_026958
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCAGCCTCCGCCATAAAGGGGCTAAGCGCACTGCGTTCCAGTACTGGTCGGCCTATTGCGTTTGTCAG
AAAAATTCTTGGACCGCGGCAGCGAGTGAGCTGAGAGAACATTTGCACAGTTTGGCCATGTAAGAAG
GTGCACTGTACCTTTTGACAAAGAGACTGGCTTTCACAGAGGCATGGGTTGGGTTTCAGTTTCTCTCAG
GAAGAAGTTCAGAAATGCACTACAACAAGAATCATATTATTGATGGAGTAAAAATCCACGTTCAAGCTC
AAGAGCAAAGGCTTTGCACGGAGCTCAAACATCTGATGAAGAAAGATTCTGAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR218296 representing NM_026958
Red=Cloning site Green=Tags(s)

MAASAIKGLSALRSSTGRPIAFVRKIPWTAASELREHFAQFGHVRRCTVPFDKETGFHRGMGWVQFSSQ
EELQNALQQEHIIIDGVKIHVQAQRAKALHGAQTSDEERFLR

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDK**V**

Species: Mouse
Serotype: AAV-2
ACCN: NM_026958
ORF Size: 336 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_026958.3</u> , <u>NP_081234.2</u>
RefSeq Size:	427 bp
RefSeq ORF:	339 bp
Locus ID:	380773
UniProt ID:	<u>Q9D8T7</u>
Cytogenetics:	12 D2
MW:	13.1 kDa