

## Product datasheet for **MR218299L4V**

### Spag6l (NM\_015773) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Spag6l (NM_015773) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Spag6l
Synonyms:	PF16; Spag6; Spag6l1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015773
ORF Size:	1524 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218299).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015773.2</a> , <a href="#">NP_056588.1</a>
RefSeq Size:	2483 bp
RefSeq ORF:	1524 bp
Locus ID:	50525
UniProt ID:	<a href="#">Q9JLI7</a>
Cytogenetics:	16 A3
Gene Summary:	Important for structural integrity of the central apparatus in the sperm tail and for flagellar motility.[UniProtKB/Swiss-Prot Function]



[View online »](#)