

Product datasheet for **MR218315L4V**

Cgref1 (NM_026770) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cgref1 (NM_026770) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cgref1
Synonyms:	1110004G24Rik; AW045201; Cgr11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026770
ORF Size:	975 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218315).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026770.1 , NP_081046.1
RefSeq Size:	1350 bp
RefSeq ORF:	978 bp
Locus ID:	68567
UniProt ID:	Q8R1U2
Cytogenetics:	5 B1
Gene Summary:	Mediates cell-cell adhesion in a calcium-dependent manner. Able to inhibit growth in several cell lines (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)