

Protein Sequence: >MR218344 representing NM_175189
 Red=Cloning site Green=Tags(s)

MKRREGALSRASRALRLSPFVYLLLIQVPVLEGVNITSPVRLIHGTVGKSALLSVQYSSTSSDKPVVKWQ
 LKRDKPVTVVQSIGTEVIGTLRPDYRDRIRLFENGSLLLSDLQLADEGTYEVEISITDDTFTGEKINLT
 VDVPISRPOVLVASTTVLELSEAFTLNCSHENGTKPSYTWLKDGPPLNDSRMLLSPDQKVLITITRVLME
 DDDL YSCVVENPISQVRS L PVKITVYRRSLYIILSTGGIFLLVTLVTVCACWKPSKKSRRKKRLEKQNS
 LEYMDQNDRLKSEADTLPRSGEQERKNPMALYILKDKDSSEPDPENPATEPRSTTEPGPPGYSVSPVPVG
 RSPGLPIRSARRYPRSPARSPATGRTHTSPPRAPSSPGRSRSSRSLRTAGVQRIREQDESGQVEISA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9018_b12.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_175189

ORF Size: 1254 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

| | |
|-------------------------------|---|
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_175189.4, NP_780398.2</u> |
| RefSeq Size: | 3103 bp |
| RefSeq ORF: | 1257 bp |
| Locus ID: | 72927 |
| UniProt ID: | <u>Q640R3</u> |
| Cytogenetics: | 9 A4 |
| MW: | 46.8 kDa |
| Gene Summary: | Involved in regulating cell motility and cell-matrix interactions. May inhibit cell growth through suppression of cell proliferation (By similarity).[UniProtKB/Swiss-Prot Function] |

Product images:


Circular map for MR218344