

## Product datasheet for **MR218346L4V**

### SytI2 (NM\_001040086) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SytI2 (NM_001040086) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	SytI2
Synonyms:	AI266830; mKIAA1597; Slp2; Slp2-a; Slp2-b; Slp2-c; Slp2-d
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001040086
ORF Size:	3840 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218346).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001040086.1</a> , <a href="#">NP_001035175.1</a>
RefSeq Size:	4843 bp
RefSeq ORF:	3843 bp
Locus ID:	83671
Cytogenetics:	7 E1



[View online »](#)

**Gene Summary:**

Isoform 11 acts as a RAB27A effector protein and plays a role in cytotoxic granule exocytosis in lymphocytes. Required for cytotoxic granule docking at the immunologic synapse. Isoform 1 may play a role in melanosome transport and vesicle trafficking. It controls melanosome distribution in the cell periphery and regulates melanocyte morphology. Isoform 1 acts as a positive mediator of mucus secretion by the surface mucus cells of the stomach. Mediates basal mucus secretion by gastric surface cells by promoting the proper granule biogenesis and docking of mucus granules with the apical plasma membrane.[UniProtKB/Swiss-Prot Function]