

## Product datasheet for **MR218401L3V**

### **Pde1a (NM\_001009979) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

<b>Product Type:</b>	Lentiviral Particles
<b>Product Name:</b>	Pde1a (NM_001009979) Mouse Tagged ORF Clone Lentiviral Particle
<b>Symbol:</b>	Pde1a
<b>Synonyms:</b>	AI987702; AW125737
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Vector:</b>	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
<b>Tag:</b>	Myc-DDK
<b>ACCN:</b>	NM_001009979
<b>ORF Size:</b>	1395 bp
<b>ORF Nucleotide Sequence:</b>	The ORF insert of this clone is exactly the same as(MR218401).
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>RefSeq:</b>	<a href="#">NM_001009979.1</a> , <a href="#">NP_001009979.1</a>
<b>RefSeq Size:</b>	1662 bp
<b>RefSeq ORF:</b>	1398 bp
<b>Locus ID:</b>	18573
<b>UniProt ID:</b>	<a href="#">Q61481</a>
<b>Cytogenetics:</b>	2 C3
<b>Gene Summary:</b>	Cyclic nucleotide phosphodiesterase with a dual-specificity for the second messengers cAMP and cGMP, which are key regulators of many important physiological processes. Has a higher affinity for cGMP than for cAMP (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)