

## Product datasheet for MR218405L3

### Atrn (NM\_009730) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Atrn (NM_009730) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Atrn
Synonyms:	AW558010; mg; Mgc; Mgca
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218405).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	NM_009730
ORF Size:	4284 bp



[View online >](#)

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_009730.2</a> , <a href="#">NP_033860.2</a>
<b>RefSeq Size:</b>	8740 bp
<b>RefSeq ORF:</b>	4287 bp
<b>Locus ID:</b>	11990
<b>UniProt ID:</b>	<a href="#">Q9WU60</a>
<b>Cytogenetics:</b>	2 63.26 cM
<b>Gene Summary:</b>	This gene encodes a widely expressed transmembrane glycoprotein that plays important roles in diverse physiological processes such as regulation of hair pigmentation, monocyte-T cell interaction and control of energy homeostasis. The encoded preproprotein undergoes proteolytic processing to generate a mature, functional protein. Certain mutations in this gene are responsible for the mahogany mouse phenotype of dark brown or black coat on a normally agouti background. Mice with loss-of-function mutations in this gene exhibit black coat color, tremor, adiposity, higher basal metabolic rate, juvenile-onset hypomyelination and age-dependent spongiform neurodegeneration of the central nervous system. [provided by RefSeq, Jul 2016]

Product images:



Circular map for MR218405L3