

## Product datasheet for **MR218405L4V**

### Atrn (NM\_009730) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Atrn (NM_009730) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Atrn
Synonyms:	AW558010; mg; Mgc; Mgca
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009730
ORF Size:	4284 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218405).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009730.2</a> , <a href="#">NP_033860.2</a>
RefSeq Size:	8740 bp
RefSeq ORF:	4287 bp
Locus ID:	11990
UniProt ID:	<a href="#">Q9WU60</a>
Cytogenetics:	2 63.26 cM



[View online »](#)

**Gene Summary:**

This gene encodes a widely expressed transmembrane glycoprotein that plays important roles in diverse physiological processes such as regulation of hair pigmentation, monocyte-T cell interaction and control of energy homeostasis. The encoded preproprotein undergoes proteolytic processing to generate a mature, functional protein. Certain mutations in this gene are responsible for the mahogany mouse phenotype of dark brown or black coat on a normally agouti background. Mice with loss-of-function mutations in this gene exhibit black coat color, tremor, adiposity, higher basal metabolic rate, juvenile-onset hypomyelination and age-dependent spongiform neurodegeneration of the central nervous system. [provided by RefSeq, Jul 2016]