



Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001033141
ORF Size:	690 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<a href="#">NM_001033141.1</a>
RefSeq Size:	1111 bp
RefSeq ORF:	693 bp
Locus ID:	68545
UniProt ID:	<a href="#">Q3TZW0</a>
Cytogenetics:	18 B2
MW:	24.9 kDa