

Product datasheet for **MR218440L4V**

Tmem132d (NM_172885) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Tmem132d
Synonyms:	C630028F04Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172885
ORF Size:	3291 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR218440).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_172885.2 , NP_766473.1
RefSeq Size:	4278 bp
RefSeq ORF:	3294 bp
Locus ID:	243274
UniProt ID:	Q76HP3
Cytogenetics:	5 G1.2-G1.3



Gene Summary:

May serve as a cell-surface marker for oligodendrocyte differentiation.[UniProtKB/Swiss-Prot Function]