

## Product datasheet for MR218444L3

### Tpm1 (NM\_001164252) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Tpm1 (NM_001164252) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Tpm1
Synonyms:	AA986836; AI854628; alpha-TM; TM2; Tm3; Tmpa; Tpm-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218444).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



**Locus ID:** 22003

**Cytogenetics:** 9 36.27 cM

**Gene Summary:** Binds to actin filaments in muscle and non-muscle cells. Plays a central role, in association with the troponin complex, in the calcium dependent regulation of vertebrate striated muscle contraction. Smooth muscle contraction is regulated by interaction with caldesmon. In non-muscle cells is implicated in stabilizing cytoskeleton actin filaments.[UniProtKB/Swiss-Prot Function]