

Product datasheet for **MR218452**

Faim2 (NM_001038658) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Faim2 (NM_001038658) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Faim2
Synonyms:	2900002L20Rik; Lfg; Lfg2; NMP25; Nmp35; Tmbim2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR218452 representing NM_001038658 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACCCAGGGAAAGCTCTCTGTGGCTAACAAGGCCCGGGACTGAGGGGCAGCAGCACCAGGCAAATG
GAGAGAAGAAGGATGCTCCGGCAGTGCCCTCGGCACCCTCCTATGAGGAGGCCACCTCTGGAGAGGG
GCTGAAGGCAGGGACCTCCCCAGGGCCCAACAGCTGTGCCACTCCATCCAAGCTGGGCCTACGTGGAC
CCCAGACATCACGAGCATTTACCACCTTCAGCTGGGATGACCAGAAAGTTCGTCGGCTCTTCATCAGAA
AGGTCTATACCATCCTATTGGTCCAGCTGCTGGTGACCTGGCTGTGGTGGCTCTTTACCTTCTGTGA
TGTGTTAAGGACTATGTCCAGGCTAACCCGGCTGGTACTGGGCGTCTATGCCGTGTTCTTTGCAACT
TACCTGACTCTGGCTTGCTGTTCTGGACCCAGGAGACATTTCCCTTGGAACCTGATTCTGCTGACCATCT
TTACCCTGTCCATGGCCTACCTCACTGGGATGCTGTCCAGTTACTACAACACCACATCTGTGCTGTGTG
CCTGGTCATCACGGCTCTCGTCTGCCTCTCGGTACCATCTTCAGTTTCCAGACCAAGTTTGACTTCACC
TCCTGTGAGGGTGTGCTGTTGTGCTGCTATGACCCCTCTTTCAGTGGACTCCTCCTCGCTGTCTCC
TGCCCTTCCAATACGTCCCCTGGCTCCACGCGGTCTATGCTGTACTGGGAGCGGGTGTGTTACGTTGTT
CCTGGCCTTTGATACCCAGTTGCTGATGGGTAACCGTCGACACTCGCTGAGCCCCGAGGAGTATATTTT
GGAGCCCTCAACATCTACTTAGACATTATCTATATCTTACCTTCTTCTGCAGCTTTTTGGCACCAACC
GGGAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR218452 representing NM_001038658
 Red=Cloning site Green=Tags(s)

MTQGKLSVANKAPGTEGQQHQANGEKKDAPVPSAPPSYEEATSGEGLKAGTFPQGPTAVPLHPSWAYVD
 PRHHEHFTTFSWDDQKVRRLFIRKVYITLLVQLLVTLAVVALFTFCDVVKDYVQANPGWYWASYAVFFAT
 YLTLACCSGPRRHFPWNLLITIFTLSMAYLTGMLSSYYNTTSVLLCLVITALVCLSVTIFSFQTKFDFT
 SCQGVLFVLLMTLFFSGLLLAVLLPFQYVPWLHAVYAVLGAGVFTLFLAFDTQLLMGNRRRHSLSPEEYIF
 GALNIYLDIIYIFTFLLQLFGTNR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001038658

ORF Size: 915 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001038658.2](#), [NP_001033747.1](#)

RefSeq Size: 4484 bp

RefSeq ORF: 918 bp

Locus ID: 72393


UniProt ID: [Q8K097](#)

Cytogenetics: 15 F1

MW: 34.7 kDa

Gene Summary: Antiapoptotic protein which protects cells uniquely from Fas-induced apoptosis. Regulates Fas-mediated apoptosis in neurons by interfering with caspase-8 activation. Plays a role in cerebellar development by affecting cerebellar size, internal granular layer (IGL) thickness, and Purkinje cell (PC) development.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR218452