

Product datasheet for **MR218480A1V**

Mouse Gng8 (NM_010320) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Gng8 (NM_010320) AAV Particle
Tag: Myc-DDK
Symbol: Gng8
Synonyms: G(y)8
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR218480 representing NM_010320
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCCAACAACATGGCCAAGATTGCTGAGGCCCGCAAGACCGTGGAGCAGCTGAAGCTGGAAGTGAACA
TCGATCGCATGAAGGTGTCGCAGGCAGCAGCCGAGCTTTGGCTTTCTGCGAAACGCACGCTAAGGATGA
CCCACTGGTGACTCCTGTCCCTGCCGCCGAGAATCCCTTCCGCGACAAGCGACTCTTTTGACCCCTGCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR218480 representing NM_010320
Red=Cloning site Green=Tags(s)

MSNNMAKIAEARKTVEQLKLEVNIDRMKVSQAAAELLAFCEHAKDDPLVTPVPAENPFRDKRLFCTLL

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_010320
ORF Size: 210 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [NM_010320.3](#), [NP_034450.1](#)

RefSeq Size: 920 bp

RefSeq ORF: 213 bp

Locus ID: 14709

UniProt ID: [P63078](#)

Cytogenetics: 7 9.15 cM

MW: 7.8 kDa