

Product datasheet for **MR218480L4V**

Gng8 (NM_010320) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gng8 (NM_010320) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gng8
Synonyms:	G(y)8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010320
ORF Size:	213 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218480).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010320.3 , NP_034450.1
RefSeq Size:	920 bp
RefSeq ORF:	213 bp
Locus ID:	14709
UniProt ID:	P63078
Cytogenetics:	7 9.15 cM



[View online »](#)

Gene Summary:

Guanine nucleotide-binding proteins (G proteins) are involved as a modulator or transducer in various transmembrane signaling systems. The beta and gamma chains are required for the GTPase activity, for replacement of GDP by GTP, and for G protein-effector interaction. This subunit may have a very specific role in the development and turnover of olfactory and vomeronasal neurons.[UniProtKB/Swiss-Prot Function]