

## Product datasheet for **MR218501L3V**

### Sfmbt1 (NM\_001166531) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sfmbt1 (NM_001166531) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sfmbt1
Synonyms:	4930442N21Rik; 9330180L21Rik; AA536974; Sfmbt; Smr
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001166531
ORF Size:	2589 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218501).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001166531.1</a> , <a href="#">NP_001160003.1</a>
RefSeq Size:	7743 bp
RefSeq ORF:	2592 bp
Locus ID:	54650
UniProt ID:	<a href="#">Q9JMD1</a>
Cytogenetics:	14 18.99 cM



[View online »](#)

**Gene Summary:**

Histone-binding protein, which is part of various corepressor complexes. Mediates the recruitment of corepressor complexes to target genes, followed by chromatin compaction and repression of transcription. Plays a role during myogenesis: required for the maintenance of undifferentiated states of myogenic progenitor cells via interaction with MYOD1. Interaction with MYOD1 leads to the recruitment of associated corepressors and silencing of MYOD1 target genes. Part of the SLC complex in germ cells, where it may play a role during spermatogenesis.[UniProtKB/Swiss-Prot Function]