

## Product datasheet for **MR218503**

### **Rft1 (NM\_177815) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Rft1 (NM_177815) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rft1
Synonyms:	D930025H04
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MR218503 representing NM\_177815  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGGGCAGCCAGGAGGTGCTAGGCCAGGCGGCCGGTTGGCTCCTCGGTCTCCTGCTGCAGTCTTGT  
 TTCGATTAATCACCTTTGCTTGAATGCATTTATTCTTCGATTCCGTCAAAGGAGATTGTGGGCATAGT  
 GAATGTCAGACTAACGCTGCTTTACTCAACCACAACCTTCTGGCCAGGAGGCCTTCCGAAGAGCTTGT  
 CTGAGCGGGGGTCCAGCGGGACTGGAGTCAGACCCTCAACCTCCTGTGGCTAACTGTCCCACTGGGTA  
 TATTTTGGTCTCCTGCTTGGGCTGGTCTGGTTGCAGTTGCTTGAAGTGCCTGATCCAGATGTGGTCCC  
 CTAATATGGAAGTGGTGTGCTGTTCTTTGGTCTCTCAGCAGTTGTGAACTTCTGGGAGAGCCTTTCTGG  
 GTCTTGGCACAAGCACATATGTTTGTAAAGCTCAAGTTCTTGTGAGAGCATGTCGGTCATCCTTAGGA  
 GTGTCTGACGGCTCTTCTGCTGTGGCTGCCTCACTGGGGTCTGTACATCTTTCTTTGGCTCAGCT  
 TCTCTACAGACTGTTCTGGTGTCTGCTATGCCATCTATTTAATAACAATTACTGCGTTCGCCAGAAATCA  
 GCCAAGCAACTAACTCTTCTGTCTCCAGAGTAACACAACCTGTTGCCAGTATTTCAAGAAGCAGAGCAT  
 TTGTCAACTGGAAGGAGGAGGATTGGCTTGGAGCTTCTCAAACAGTCGTTCTTGAAGCAGATTTTGAC  
 AGAAGGCGAGCGGTACGTGATGACATTTTTGAATGTGTTAAACTTTGGCGATCAGGGTGTATATGACATA  
 GTGAATAATCTTGGCTCCCTTGTGGCCAGGCTAATTTTCCAGCCAGTAGAGGAGAGTTTTTATTTATTCT  
 TTGCTAAAGTGTAGAGAGAGAGAAGGACGCCTCACTTCAAGAAGCAGGACGATGTTGCTGTGGCTGCTGC  
 GGTTTTGGAGTCCCTCCTCAAGCTGGCCTTGTGACTGGACTGACCATGACTGTTTTTGGCTTTGCCTAT  
 TCCCAACTGGCTTTGGACATCTATGGAGGAGCCATGCTTAGTTCGGGACTGCGCCCTGTGCTGATGCGTT  
 GCTATTGCCTCTATGCTCCTGCTTGTATCAATGGTGTGACCGAGTGTTCATGTTTGGCCCATGAG  
 CAAAGAGGAGGTTGACAGGTACAATTTACCATGCTGGCACTGTATCCTCGTTTTTAGTCTTATCCTAC  
 CTGCTGACCAGTTGGTGTGGCAGTGTGGCTTTCATCATGGCCAACGCTTTAACATGGGCATCCGGATCA  
 CACAGAGTCTTTCCTTCCATCACTACTTCCGAGAGAGCCCTCACAGACCCCTGGCTGGCCTGCGCCT  
 GTCACCAGTCTTCTTGGGGTGTTCATTCTCAGTGTGGGATTACTAGTGTCTCAGAGGCATTCTCTGCT  
 TGTGAGAGGGGCTGGCCAGCTAGGCTGGCACACATCGTGTGGGGACCATCTGCTTAGGAGTGACTCTTG  
 GGACAGCGTTCCTCACAGAGACCAAGCTGATCCACTTTCTCAGGACTCAGTTGGGTAGATCCAGACTCAG  
 TGACAAAATGACA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**

>MR218503 representing NM\_177815  
 Red=Cloning site Green=Tags(s)

MGSQEVLGQAARLASSGLLLQVLFRLITFVNAFILRFLSKEIVGIVNVRLTLLYSTTTFLAREAFRRAC  
 LSGGAQRDWSQTLNLLWLVPLGIFWSSCLGWVWLQLEVPDPDVPYVYGTGVLFFGLSAVVELLGEPFW  
 VLAQAHMFVKLVLAESMSVILRSVLTALLVLWLPHWGLYIFSLAQLLYTTVLVLCYAIYLIQLLRSPES  
 AKQLTLPVSRVTQLLPSISRSRAFVNWKEAGLAWSFFKQSFLKQILTEGERYVMTFLNVLNFGDQGVYDI  
 VNNLGSVARLIFQPVEESFYLFFAKVLEREKDASLQKQDDVAVAAAVLESLLKALLTGLTMTVFGFAY  
 SQLALDIYGGAMLSSGSGPVLRCYCLYVLLLAINGVTECFMFAAMSKEEVDRYNFMLALSSSFLVLSY  
 LLTSWCGSVGFIMANCFNMGIRITQSLSFIIHYFRESPHRPLAGLRSPVLLGVFVILSAGITSVSEAFLC  
 CERGWPARLAHIAVGTICLVTLGTAFLTETKLIHFLRTQLGRSRLSDKMT

**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV

**Chromatograms:**

[https://cdn.origene.com/chromatograms/mm9015\\_e11.zip](https://cdn.origene.com/chromatograms/mm9015_e11.zip)

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**

**ACCN:** NM\_177815

**ORF Size:** 1623 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_177815.3](#), [NP\\_808483.2](#)
**RefSeq Size:** 2384 bp

**RefSeq ORF:** 1626 bp

**Locus ID:** 328370


**UniProt ID:** [Q8C3B8](#)

**Cytogenetics:** 14 B

**MW:** 60.8 kDa

**Gene Summary:** May be involved in N-linked oligosaccharide assembly. May participate in the translocation of oligosaccharide from the cytoplasmic side to the luminal side of the endoplasmic reticulum membrane (By similarity).[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR218503