

Product datasheet for **MR218523L3V**

Lactb (NM_030717) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lactb (NM_030717) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lactb
Synonyms:	LACT-1; Lact1; Mrpl56
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_030717
ORF Size:	1653 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218523).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030717.1 , NP_109642.1
RefSeq Size:	2031 bp
RefSeq ORF:	1656 bp
Locus ID:	80907
UniProt ID:	Q9EP89
Cytogenetics:	9 C



[View online »](#)

Gene Summary:

Mitochondrial serine protease that acts as a regulator of mitochondrial lipid metabolism (By similarity). Acts by decreasing protein levels of PISD, a mitochondrial enzyme that converts phosphatidylserine (PtdSer) to phosphatidylethanolamine (PtdEtn), thereby affecting mitochondrial lipid metabolism (By similarity). It is unclear whether it acts directly by mediating proteolysis of PISD or by mediating proteolysis of another lipid metabolism protein (By similarity). Acts as a tumor suppressor that has the ability to inhibit proliferation of multiple types of cancer cells: probably by promoting decreased levels of PISD, thereby affecting mitochondrial lipid metabolism (PubMed:28329758).[UniProtKB/Swiss-Prot Function]