

Product datasheet for **MR218524L3V**

Ctnna2 (NM_009819) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ctnna2 (NM_009819) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ctnna2
Synonyms:	AI481747; Catna; Catna2; cdf; chp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009819
ORF Size:	2859 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218524).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009819.2 , NP_033949.2
RefSeq Size:	4046 bp
RefSeq ORF:	2862 bp
Locus ID:	12386
UniProt ID:	Q61301
Cytogenetics:	6 33.54 cM



[View online »](#)

Gene Summary:

May function as a linker between cadherin adhesion receptors and the cytoskeleton to regulate cell-cell adhesion and differentiation in the nervous system (PubMed:12123610, PubMed:15034585). Required for proper regulation of cortical neuronal migration and neurite growth. It acts as negative regulator of Arp2/3 complex activity and Arp2/3-mediated actin polymerization. It thereby suppresses excessive actin branching which would impair neurite growth and stability (By similarity). Regulates morphological plasticity of synapses and cerebellar and hippocampal lamination during development. Functions in the control of startle modulation (PubMed:12089526).[UniProtKB/Swiss-Prot Function]