

## Product datasheet for **MR218534L4V**

### **Nucks1 (NM\_001145804) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Nucks1 (NM_001145804) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nucks1
Synonyms:	2700010L10Rik; 8430423A01Rik; A1647518; C78391; Nucks
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001145804
ORF Size:	699 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218534).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001145804.1</a> , <a href="#">NP_001139276.1</a>
RefSeq Size:	6093 bp
RefSeq ORF:	702 bp
Locus ID:	98415
UniProt ID:	<a href="#">Q80XU3</a>
Cytogenetics:	1 E4



[View online »](#)

**Gene Summary:**

Chromatin-associated protein involved in DNA repair by promoting homologous recombination (HR). Binds double-stranded DNA (dsDNA) and secondary DNA structures, such as D-loop structures, but with less affinity than RAD51AP1.[UniProtKB/Swiss-Prot Function]