

Product datasheet for MR218537

Klrb1a (NM_001159902) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Klrb1a (NM_001159902) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Klrb1a
Synonyms: ly-55A; Ly55a; NKR-P1.7; NKR-P1A; Nkrp1-a; Nkrp1a; NKR12
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR218537 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGCATCTCTATGCACAATGGACACAGCAAGGGTCTACTTTGGTTAAAGCCACCCAGGACTCCAGGGG
CTTGGCATGAGTACCCCATCTCTCCCCAGTGGAGTCTAATAACAAAACCATCAATAGAAAATG
CTATGTGCTTATTCAAGAGAACCTGAATAAAACAACAGATTGTTCAAGCTAAGCTAGAGTGCCACAAGAC
TGGCTTTCACACCGAGATAAGTGCTTTCACGTTTCTCAAGTTTCCAACACTTGGGAGGAAGTCTAGTTG
ATTGTGATGGAAAAGGAGCCACTTTGATGCTCATTCAAGACCAAGAAGAACTGAGATTCCTACTGGACTC
AATAAAGGAAAAATACAATTCATTTGGATTGGACTAAGGTACACATTGCCAGACATGAAGTGAAGTGG
ATAAATGGATCGACTTTAAATTCGATGTATTAATAACTACTGGTACACTGAAAATGACAGCTGTGCTG
CTATCTCAGGAGACAAAGTACTTTGAGAGCTGTAATTCAGACAACCGTTGGATCTGCCAAAAGGAACT
ATACCATGAAACCTGAGCAACTATGTGGGTTATGGACAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR218537 protein sequence
Red=Cloning site Green=Tags(s)

MHLLCTMDTARVYFGLKPPRTPGAWHESPPSLPPVRVLIQKPSIEKCYVLIQENLNKTTDCSAKLECPQD
WLSHRDKCFHVSQVSNTWEEGLVDCDGGATLMLIQDQEELRFLLDSEIKEYNSFWIGLRYTLPDMNWKW
INGSTLNSDVLKITGDTENDSCAAISGDKVTFESCNSDNRWICQKELYHETLSNYVGYGH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001159902

ORF Size: 600 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq Size: 1185 bp

RefSeq ORF: 603 bp

Locus ID: 17057

UniProt ID: [P27811](#)

Cytogenetics: 6 63.06 cM

MW: 22.9 kDa

Gene Summary: Plays a stimulatory role on natural killer (NK) cell cytotoxicity.[UniProtKB/Swiss-Prot Function]