

Product datasheet for **MR218555L3V**

C5ar2 (NM_001146005) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C5ar2 (NM_001146005) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	C5ar2
Synonyms:	C5L2; E030029A11Rik; Gpr77
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001146005
ORF Size:	1074 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218555).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001146005.1 , NP_001139477.1
RefSeq Size:	3429 bp
RefSeq ORF:	1077 bp
Locus ID:	319430
UniProt ID:	Q8BW93
Cytogenetics:	7 A2



[View online »](#)

Gene Summary:

Receptor for the chemotactic and inflammatory C3a, C4a and C5a anaphylatoxin peptides and also for their dearginated forms ASP/C3adesArg, C4adesArg and C5adesArg respectively. Couples weakly to G(i)-mediated signaling pathways (By similarity).[UniProtKB/Swiss-Prot Function]