

Product datasheet for **MR218556**

Lin7a (NM_001039354) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lin7a (NM_001039354) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Lin7a
Synonyms:	A1848705; LIN-7A; MALS-1; TIP-33; Veli; Veli1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR218556 representing NM_001039354 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGAAGCCGAGCGTCACTTCGGCTCCACGGCAGACATGGCGACATTGACAGTGGTCCAGCCGCTTA
CTCTGGACAGAGATGTTGCAAGCAATCGAACTACTGGAAAAGCTACAAGAATCTGGAGAAGTACCAGT
GCACAAGCTCCAGTCTCTCAAAAAGGTGCTTCAGAGTGAGTTTTGTACAGCAATCCGAGAGGTGTATCAA
TACATGCATGAAACGATTACTGTTAATGGCTGCCCTGAATTCGTGCGAGGGCCACAGCAAAGGCAACAG
TTGCAGCTTTTGCAGCCAGCGAAGGCCACTCCACCCCTCGGGTAGTTGAACTGCCAAAGACTGATGAAGG
CCTGGGTTTTAATGTGATGGGAGGGAAGGAACAGAATCCCGATTTACATCTCCCGCATCATCCCTGGA
GGGGTGGCTGAAAGACATGGAGGCTCAAAGAGGGGACCAAGCTGCTATCAGTGAACGGAGTGAGTGTGG
AAGGGGAGCACCATGAGAAAGCTGTGAACTTCTCAAGGCTGCTAAGGACAGTGTGAAGCTGGTGGTCA
ATACACCCAAAAGTCTGGAAGAGATGGAGGCTCGTTTCGAGAAGCTGCGGACAGCTCGGCGTCGACAG
CAGCAGCAGTTGCTATTCAGCAGCAGCAACAGCAGCAGCAACAACAACCAACAAAACCAATGTCA

ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR218556 representing NM_001039354
Red=Cloning site Green=Tags(s)

MLKPSVTSAPTADMATLTVVQPLTLDRDVARAIELLEKLQESGEVPVHKLQSLKVKVQSEFCTAIREVYQ
 YMHETITVNGCPEFRARATAKATVAFAASEGHSHPRVVELPKTDEGLGFNVMGGKEQNSPIYISRIIPG
 GVAERHGGLKRGDQLLSVNGVSVEGEHHEKAVELLKAAKDSVKLVVRYTPKVLEEMEARFEKLRTARRRQ
 QQQLLIQQQQQQQQPQQNHMS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9080_d12.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001039354

ORF Size: 699 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001039354.1](#), [NP_001034443.1](#)

RefSeq Size: 5887 bp

RefSeq ORF: 702 bp

Locus ID: 108030


UniProt ID: [Q8JZS0](#)

Cytogenetics: 10 D1

MW: 26 kDa

Gene Summary: Plays a role in establishing and maintaining the asymmetric distribution of channels and receptors at the plasma membrane of polarized cells. Forms membrane-associated multiprotein complexes that may regulate delivery and recycling of proteins to the correct membrane domains. The tripartite complex composed of LIN7 (LIN7A, LIN7B or LIN7C), CASK and APBA1 may have the potential to couple synaptic vesicle exocytosis to cell adhesion in brain. Ensures the proper localization of GRIN2B (subunit 2B of the NMDA receptor) to neuronal postsynaptic density and may function in localizing synaptic vesicles at synapses where it is recruited by beta-catenin and cadherin. Required to localize Kir2 channels, GABA transporter (SLC6A12) and EGFR/ERBB1, ERBB2, ERBB3 and ERBB4 to the basolateral membrane of epithelial cells.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR218556