

## Product datasheet for **MR218583L4V**

### Ramp1 (NM\_178401) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ramp1 (NM_178401) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ramp1
Synonyms:	9130218E19Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178401
ORF Size:	309 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218583).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_178401.3</a> , <a href="#">NP_848488.2</a>
RefSeq Size:	868 bp
RefSeq ORF:	312 bp
Locus ID:	51801
Cytogenetics:	1 46.08 cM
Gene Summary:	Transports the calcitonin gene-related peptide type 1 receptor (CALCRL) to the plasma membrane. Acts as a receptor for calcitonin-gene-related peptide (CGRP) together with CALCRL.[UniProtKB/Swiss-Prot Function]



[View online »](#)