

Product datasheet for **MR218594**

Fxyd4 (NM_001173372) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fxyd4 (NM_001173372) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Fxyd4
Synonyms: 0610008I02Rik; A1267073; Chi; Chif
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR218594 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGGAAATCACCTGTGCCTTCTCCTGCTGCTAGCAGGTCTGCCGGCCTTGAAGCCAGTGACCCAG
 TTGATAAAGACAGTCCCTTCTACTATGACTGGGAGAGCCTGCAGCTGGGAGGATTGATTTTTGGAGGGCT
 CCTGTGCATCGTGGAATTGCCATGGCCCTGAGTGGAAGTGCAAATGTAGGCGTACCCATAAGCCCAGT
 TCCTTACCTGAAAAGCCACTCCACTCATCTCCAGGCTCTGCCAATACCTGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >MR218594 protein sequence
 Red=Cloning site Green=Tags(s)
 MEEITCAFLLLL LAGLPAL EASDPVDK DSPFYDWESLQLGGLIFGLLCIAGIAMALSGKCKCRRT HKPS
 SLP GKATPLIIPGSANTC

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001173372

ORF Size: 267 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001173372.1](#), [NP_001166843.1](#)

RefSeq Size: 551 bp

RefSeq ORF: 267 bp

Locus ID: 108017


UniProt ID: [Q9D2W0](#)

Cytogenetics: 6 F1

MW: 9.3 kDa

Gene Summary: This gene encodes a member of a family of small membrane proteins that share a 35-amino acid signature sequence domain, beginning with the sequence PFXYD and containing 7 invariant and 6 highly conserved amino acids. The approved human gene nomenclature for the family is FXYD-domain containing ion transport regulator. Mouse FXYD5 has been termed RIC (Related to Ion Channel). FXYD2, also known as the gamma subunit of the Na,K-ATPase, regulates the properties of that enzyme. FXYD1 (phospholemman), FXYD2 (gamma), FXYD3 (MAT-8), FXYD4 (CHIF), and FXYD5 (RIC) have been shown to induce channel activity in experimental expression systems. Transmembrane topology has been established for two family members (FXYD1 and FXYD2), with the N-terminus extracellular and the C-terminus on the cytoplasmic side of the membrane. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR218594