

## Product datasheet for **MR218595L3V**

### Pigx (NM\_001111025) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Pigx (NM_001111025) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pigx
Synonyms:	2010319C14Rik; PIG-X
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001111025
ORF Size:	756 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218595).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001111025.1</a> , <a href="#">NP_001104495.1</a>
RefSeq Size:	1013 bp
RefSeq ORF:	759 bp
Locus ID:	72084
Cytogenetics:	16 B2



[View online »](#)

**Gene Summary:**

This gene encodes a type I transmembrane protein in the endoplasmic reticulum (ER). The protein is an essential component of glycosylphosphatidylinositol-mannosyltransferase I, which transfers the first of the four mannoses in the GPI-anchor precursors during GPI-anchor biosynthesis. Studies in rat indicate that the protein is translated from a non-AUG translation initiation site. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]