

Product datasheet for **MR218598L4V**

Mrpl47 (NM_029017) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mrpl47 (NM_029017) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mrpl47
Synonyms:	4833424P18Rik; CGI-20; CGI-204; Gm9859; L47mt; MRP-L47; MTF/L47; NCM; NCM1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_029017
ORF Size:	759 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218598).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029017.2 , NP_083293.1
RefSeq Size:	899 bp
RefSeq ORF:	759 bp
Locus ID:	74600
UniProt ID:	Q8K2Y7
Cytogenetics:	3 A3


[View online »](#)

Gene Summary:

Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. This gene is immediately adjacent to the gene for BRG1/brm-associated factor 53A (also known as BAF complex 53 kDa subunit protein A in humans) in a tail-to-tail orientation. [provided by RefSeq, Jul 2008]