

Product datasheet for **MR218610L4V**

Veph1 (NM_145820) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Veph1 (NM_145820) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Veph1
Synonyms:	2810471M23Rik; Veph
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145820
ORF Size:	2499 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218610).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145820.3 , NP_665819.2
RefSeq Size:	6450 bp
RefSeq ORF:	2502 bp
Locus ID:	72789
UniProt ID:	A1A535
Cytogenetics:	3 E1
Gene Summary:	Interacts with TGF-beta receptor type-1 (TGFBR1) and inhibits dissociation of activated SMAD2 from TGFBR1, impeding its nuclear accumulation and resulting in impaired TGF-beta signaling. May also affect FOXO, Hippo and Wnt signaling.[UniProtKB/Swiss-Prot Function]



[View online »](#)