

## Product datasheet for **MR218613L4V**

### Galnt6 (NM\_001161768) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Galnt6 (NM_001161768) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Galnt6
Synonyms:	4632410F13; AW047994; GalNAc-T6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001161768
ORF Size:	1869 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218613).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001161768.1</a> , <a href="#">NP_001155240.1</a>
RefSeq Size:	5712 bp
RefSeq ORF:	1869 bp
Locus ID:	207839
UniProt ID:	<a href="#">Q8C7U7</a>
Cytogenetics:	15 F1



[View online »](#)

**Gene Summary:**

Catalyzes the initial reaction in O-linked oligosaccharide biosynthesis, the transfer of an N-acetyl-D-galactosamine residue to a serine or threonine residue on the protein receptor. May participate in synthesis of oncofetal fibronectin. Has activity toward Muc1a, Muc2, EA2 and fibronectin peptides (By similarity).[UniProtKB/Swiss-Prot Function]