

Product datasheet for **MR218615L4V**

Nek3 (NM_001162947) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nek3 (NM_001162947) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nek3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001162947
ORF Size:	1533 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR218615).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001162947.1 , NP_001156419.1
RefSeq Size:	1989 bp
RefSeq ORF:	1536 bp
Locus ID:	23954
UniProt ID:	Q9R0A5
Cytogenetics:	8 11.07 cM
Gene Summary:	Protein kinase which influences neuronal morphogenesis and polarity through effects on microtubules. Regulates microtubule acetylation in neurons. Contributes to prolactin-mediated phosphorylation of PXN and VAV2.[UniProtKB/Swiss-Prot Function]



[View online »](#)